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
Dem  
Nano  
Regio  
Leben zuliebe.

## INSTRUCTION FOR USE


**INTRODUCTION**  
In December 2019, the novel respiratory disease (COVID-19) caused by the coronavirus (SARS-CoV-2) was reported in Wuhan, China. According to WHO, most of the people infected with SARS-CoV-2 have mild to moderate respiratory symptoms. The detection of these antigens indicates a high probability of infection with the novel coronavirus.

**POSITIVE TEST RESULT**  
If a colored line is visible in the Test-line area (T) and a colored line in the Control line (C) the test result is positive. The test result means that SARS-CoV-2 antigen are detectable in your saliva sample. The detection of these antigens indicates with a high probability of infection with the novel coronavirus.

**Please stay at home and contact your physician or responsible authority on how to proceed further.**  
Note: The thickness of the line is insignificant, any reddish color in the Test line (T) should be considered a positive result. The positive test result must be confirmed by PCR.

**NEGATIVE TEST RESULT**  
If only a colored line is visible in the Control line area (C) the test result is negative.

**INVALID TEST RESULT**  
If there is no colored line in the Control line (C) and no colored line in the Test line (T) the test result is invalid. This indicates that the test is not working properly. The test result must be confirmed by the physician and get information on how to proceed further.

An explanation of how to read and interpret the test results requires a thorough reading of the instructions for use. It is important to fully understand the entire instruction pamphlet before proceeding with the test.

**TEST CONTENTS**  
• 1 test cassette with 1 pipette in a sealed pouch (SARS-CoV-2 Ag)  
• 1 sample tube with pre-filled sample extraction buffer  
• 1 saliva collector  
• 1 instruction pamphlet  
• 1 timer

**ADDITIONALLY REQUIRED:**  
• Clean, dry and flat surface.

**TEST PREPARATION**  
Let test cassette and components stand at a room temperature (15-30°C) before mixing the test. Lay all the supplied materials on a clean, dry and flat surface.

**TEST PERFORMANCE**  
Read the instructions for use completely before performing the test. A step-by-step instruction is given on the next page and describes the test procedure.

**EVALUATION OF TEST RESULTS**  
To read the test results simply determine whether a line is present or absent at the Control (C) position. It does not matter how strong or weak a Control line (C) is.

## STEP-BY-STEP INSTRUCTION

### STEP 1

Do not eat, drink, smoke or chew gum for 30 minutes before collecting saliva sample. This may affect the test results.


Attention: This tube does not contain liquid!

### STEP 2


Open the sealed pouch and remove the test cassette. Lay it face up on a clean, dry and flat surface.

### STEP 3

Screw saliva collection funnel on to the saliva collector tube.

### STEP 4

Place the tip of tongue against the mucous membrane of the oral cavity to emit saliva, and then swallow saliva to the throat without ingesting saliva to esophagus, followed by returning saliva to the mouth. Repeat the process for 5 times, then spit saliva into a saliva collector. Avoid bubbles while collecting saliva.

### STEP 5

Vertically hold the saliva collector tube and unscrew the saliva collector funnel to remove it from the tube. Discard the funnel.

Take sample tube with pre-filled sample extraction buffer and remove the white cap of the sample tube.

### STEP 6

Draw 2/3 dropper range of saliva sample from the saliva collector tube with the help of pipette and transfer the entire contents into the extraction buffer of the sample tube.

To thoroughly mix the saliva sample with extraction buffer, draw the solution up and down with the pipette about 10 times!

Attention: Please avoid bubbles!

### STEP 7

Close the cap of the sample tube.

### STEP 8


Add 3 full drops of the mixed solution vertically into the sample well (S) of the test cassette.

Please note that there should be no liquid applied to the result windows marked with the letters (C) and (T) or to the sample well (S) after adding the drops to the sample well (S).


Read the result 15-20 minutes after adding the sample. Result got after 20 minutes is invalid.

## For self-testing

### RESULTS

**POSITIVE**  
SARS-CoV-2 Ag

**NEGATIVE**  
SARS-CoV-2 Ag

**INVALID**  
SARS-CoV-2 Ag

**POSITIVE**  
SARS-CoV-2 Ag

**NEGATIVE**  
SARS-CoV-2 Ag

## PERFORMANCE EVALUATION

Zuhause TEST Corona	Reference tests PCR	
	Positive	Negative
Positive	99	4
Negative	6	396
Total	105	400
	505	505

Specificity: 94.29 %  
Sensitivity: 99.00 %  
Precision: 98.02 %  
The limit of detection of Zuhause test Corona has been evaluated at 1.02 x 10<sup>4</sup> TCID<sub>50</sub>/ml.

No cross reactivity with the human coronaviruses including OC-43, NL63, SARS and MERS, influenza A & B viruses, rhinoviruses and other common viruses and bacteria was observed when tested using the Zuhause TEST Corona.

### LIMIT OF DETECTION AND LITERATURE

1.) Nanshan Chen\*, Min Zhou\*, Xuan Dong\*, Jieming Qu\*, Fengyun Gong, Yang Han, Yang Qiu, Jing Wang, Ying Liu, Juan Wei, Jitao Xia, Zhongming Chen, et al. (2020) Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. LANCET. January 29, 2020.

2.) World Health Organization (Coronavirus disease 2019) [https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-\(covid-2019\)-and-the-virus-that-causes-it](https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-(covid-2019)-and-the-virus-that-causes-it) (Zugriff am 27.03.2020)

3.) World Health Organization (Coronavirus disease 2019) <https://www.who.int/health-topics/coronavirus/nCoV-19> (Zugriff am 27.03.2020)

4.) Saliva as a Candidate for COVID-19 Diagnostic Testing: A Meta-Analysis. Laszlo Mark Czumbel, Szabolcs Kiss, Melli Forjós, Ivan Doherty, Zoltán Székely, Zoltán Székely, Péter Hegyi, Martin C. Stewart and Gábor Varga

• Not following the exact instructions can affect the outcome of the test. The final diagnosis must be confirmed by a physician.

• Do not use the test if the packaging is damaged. Do not use broken test components.

• All test components are only intended to be used for this test. Do not reuse the test or test components.

• The test should be carried out immediately or within one hour after opening the foil pouch (20-30°C, humidity <60%).

• After human saliva sample is collected, the test should be used immediately or within 1 hour. If it cannot be tested immediately, it can be stored at 2-8°C for 72 hours.

• Poor vision, color blindness or poor lighting may affect your ability to interpret the test correctly.

• Refer to the local regulation in force regarding the disposal of the test components.

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	Follow instructions		In vitro diagnostic medical device (for external use)		Best before (See imprint on package)
	Store at 2-8°C, 30°C Do not freeze		Contains sufficient for 1 test		Do not reuse
	Manufacturer		Keep dry		Batch number (See imprint on package)
	Reference number				

[www.zuhausestest.de](http://www.zuhausestest.de)

REF: B60200

**CE 0483**

Instructions English  
Revision from 2021-02 (Rev. 01)

NanoPro AG  
Untergasse 8  
D-35357 Marburg  
[www.nanoopro.com](http://www.nanoopro.com)

For the  
love of life.